



# WED-2000

## Digital Ultrasound Scanner



Palmtop WED-2000, a giant in miniature, which weighs only 600g along with the probe and internal rechargeable battery, is the smallest and lightest ultrasound scanner in the world. Its sharp image, heavy duty rechargeable battery, unique shape and low weight make it competent in the emergency treatment applications.



# WED-2000

A GIANT IN MINIATURE

## Digital Ultrasound Scanner



### System Features:

5.0 TFT LCD

Net weight:600g

Dimension:148mm×112mm×42mm

Low power consumption

Battery working time for 3h continuously

Easy-to-use GUI menu

Real time/freeze/zoom display modes



### Waterproof Probe:



### Probe with FREEZE Key:



### Rechargeable Battery and Car Charger:



### Technical Specification:

Monitor: 5.0inch TFT LCD

Weight: 600g

Probes: 2.5/3.5/5.0MHz multi-frequency mechanical probe

Battery: Rechargeable/3h continuously

Image Mode: B,B+B,B+M,M

Image Magnification: 8 levels zoom in

Image Permanent Storage: 128 frames

Gray Shades: 256 levels

Character Display: Hospital, name, age, sex

Measurement: Distance, circumference, area, volume, heart rate, GA,FW,EDD

Standard Accessories: One probe

Internal rechargeable battery

Mouse

Sp-2 AC Adapter

Optional Accessories: Car charger

Sp-3 AC Adapter

Extra battery with charging cable

### Sp-3 AC Adapter:



The scanner can also be used with AC supply when it is connected to Sp-3 AC Adapter.



### Application:

We can provide specific applications according to customer's requirements.

WED-2000 is applicable to human medicine.

WED-2000V is applicable to veterinary medicine.



Shenzhen WELL.D Medical Electronics Co.,Ltd.

<http://www.welld.net> <http://www.welld.com.cn> Email: [export@welld.com.cn](mailto:export@welld.com.cn)

ADD:12/F BLock North, New Energy Bldg., Nanhai Ave., Nanshan District, ShenZhen 518054,Guangdong, China

Tel:+86-755-26072063/805/821 Fax:+86-755-26073350

Copyright (c) 2009 SZWELLD. All right reserved.

Version:2000-V3.0